

# Flexural strength of veneering ceramics for zirconia

# J. Fischer\*, B. Stawarczyk, C.H.F. Hämmerle

Clinic for Fixed and Removable Prosthodontics, Center for Dental and Oral Medicine, University of Zurich, Plattenstrasse 11, CH-8032 Zurich, Switzerland

#### ARTICLE INFO

Article history: Received 21 September 2007 Received in revised form 18 January 2008 Accepted 25 January 2008

Keywords: Zirconia Veneering ceramics Three-point flexural strength Four-point flexural strength Biaxial flexural strength

# ABSTRACT

*Objectives*: The flexural strengths of veneering ceramics for zirconia were compared. *Methods*: With 10 different veneering ceramics for zirconia (test group) and three different veneering ceramics for the metal-ceramic technique (control group) three-point flexural strength and biaxial flexural strength according to ISO 6872: 1995 as well as four-point flexural strength according to EN 843-1: 2005 were measured (n = 10). Statistical analysis was performed with one-way ANOVA and post hoc Scheffé test (SPSS, p < 0.05). *Results*: For the test group, three-point flexural strength ranged between 77.8  $\pm$  8.7 and

Results. For the test group, three-point flexural strength ranged between 77.5  $\pm$  8.7 and 106.6  $\pm$  12.5 MPa without any statistically significant differences, biaxial flexural strength between 69.1  $\pm$  4.8 and 101.4  $\pm$  10.5 MPa with three homogeneous groups and four-point flexural strength between 59.5  $\pm$  6.2 and 89.2  $\pm$  9.5 MPa with five homogeneous groups. The control group showed three-point flexural strength values ranging from 93.3  $\pm$  13.5 to 149.4  $\pm$  20.5 MPa, biaxial flexural strength values from 93.4  $\pm$  10.0 to 141.2  $\pm$  11.6 MPa, and four-point flexural strength values from 82.7  $\pm$  8.5 to 116.9  $\pm$  9.8 MPa. In every case, the results of the four-point flexure test were significantly lower than those obtained in the three-point flexure test. The three-point flexural strengths of the test group are similar to those of two ceramics of the control group. The flexural strength of one ceramic of the control group significantly exceeded the strengths of all other ceramics investigated. Conclusion: Three-point flexural strength values of veneering ceramics for zirconia are

similar to those of veneering ceramics for the metal-ceramic technique. The four-point flexure test among all three tests showed highest discrimination between the different ceramic materials.

© 2008 Elsevier Ltd. All rights reserved.

# 1. Introduction

Yttria-stabilized zirconia (Y-TZP) provides a sufficient mechanical strength to be used in frameworks for all-ceramic fixed partial dentures.<sup>1,2</sup> For esthetical reasons, these frameworks have to be veneered with an appropriate veneering ceramic. In clinical application, the veneering ceramic revealed to be the weakest link in such reconstructions.<sup>3–5</sup> Chipping of the veneer is described to be the most frequent reason for failure with a failure rate of 15.2% after a service time of  $35.1 \pm 13.8$  months.<sup>5</sup>

Among other reasons failure of a veneer may be caused by insufficient bond strength,<sup>6–8</sup> excessive tensile stress due to a thermal mismatch between veneer and framework<sup>9</sup> or excessive load due to premature contacts.<sup>10</sup> The bond strength was intensely investigated.<sup>9,11–13</sup> It revealed to be in the range of that measured with metal-ceramic systems. The tensile stress in the veneering ceramic is established during cooling after firing, when an unequal thermal contraction of both layers happens. The coefficients of thermal expansion should be adjusted in a way that during cooling a slight compression of the veneering ceramic occurs to enhance its strength.<sup>14</sup> In

<sup>\*</sup> Corresponding author. Tel.: +41 44 634 33 67; fax: +41 44 634 43 05. E-mail address: jens.fischer@zzmk.uzh.ch (J. Fischer).

<sup>0300-5712/\$ –</sup> see front matter © 2008 Elsevier Ltd. All rights reserved. doi:10.1016/j.jdent.2008.01.017

metal-ceramic systems, excessive stress to some extent may be compensated by thermal creep of the alloy, i.e. plastic flow, especially if a high gold alloy is used.<sup>15,16</sup> In all-ceramic systems, the ceramic framework is rigid and does not yield to the stress induced by a thermal mismatch to that extent. Therefore, the risk of destructive stress formed in the veneer layer might be higher in all-ceramic systems and thus would require a high mechanical strength for veneering materials for all-ceramic systems. Hence, the strength of the veneering ceramic is a crucial parameter for the clinical long-term success. For metal-ceramic restorations failure rates after 5 years, caused by chipping of the veneer are reported to be 0.4% for single  $\rm crowns^{17}$  and 2.9% for fixed partial dentures.  $^{18}$ Hence, veneering ceramics for zirconia should at least show a flexural strength, which is similar to that of veneering ceramics for alloys.

Flexural strength can be measured in a three-point flexure test, a four-point flexure test or a biaxial flexure test. In all cases, static load is applied until failure. In the three-point flexure test, a non-uniform central stress field is created, while in the four-point flexure test the stress field is uniform between the two loading pistons. In the biaxial flexure test, where a disk is loaded in the center, the probability of edge failures is reduced.<sup>19</sup> The results of the three-point flexure test and the four-point flexure test are correlated.<sup>20</sup> Lower values were found for the four-point flexure test compared to both other tests, but the relation between three-point flexure test and biaxial flexure test was not uniform for all ceramics investigated.

To the knowledge of the investigators, no systematic investigation of the flexural strength of veneering ceramics for zirconia is available.

Aim of the present study therefore was to measure the flexural strength of a variety of commercially available

veneering ceramics for zirconia to provide a comprehensive analysis of the mechanical strength of these products.

# 2. Materials and methods

Three-point flexural strength, four-point flexural strength and biaxial flexural strength of 10 different veneering ceramics for zirconia according to Table 1 were measured. As control three ceramics for the metal-ceramic technique were additionally included (Imagine Reflex, IPS d.sign, and VM13).

Specimens were prepared according to ISO 6872: 1995 (threepoint and biaxial flexural strength) or DIN EN 843-1: 2005 (fourpoint flexural strength). Separable steel molds were used to layer the ceramic. Ceramic powder and an appropriate amount of the respective liquid were mixed to form a sticky slurry, which was filled into the mold. Excess liquid was sucked off with a tissue. Only dentin was layered. Firing of the specimens was performed in a ceramic oven (Austromat D4, Dekema, Freilassing, Germany) according to the recommendations of the manufacturers (Table 2). The specimens were placed on a tray, which was covered with a layer of silica powder. After firing, the specimens were ground to the final dimensions using SiC discs P220, P500 and P1200 according to ISO 6344-1: 1998. As required by the standards the two faces of the specimens did not differ more than 0.05 mm in parallelism. Ten specimens were prepared for each series. The dimensions of the samples were measured to the next 0.01 mm. The specimens were placed in the appropriate sample holder and loaded in a universal testing machine (Z010, Zwick, Ulm, Germany) with a cross-head speed of 1 mm/min until failure. The flexural strength was calculated as mean of the 10 results.

Statistical analysis between different test methods and between the ceramics were analyzed with one-way ANOVA,

| Table 1 – Veneering ceramics used in the investigation |  |  |  |  |  |  |
|--|--|--|--|--|--|--|
|  |  |  |  |  |  |  |
| Veneering Ceramic                                      | Manufacturer                               |  |  |  |  |  |
| Cerabien ZR  | Noritake, Nagoya, Japan                    |  |  |  |  |  |
| Creation ZI  | Metalordental, Oensingen, Switzerland      |  |  |  |  |  |
| IPS e.max  | Ivoclar-Vivadent, Schaan, Liechtenstein    |  |  |  |  |  |
| Initial ZR   | GC, Tokyo, Japan                           |  |  |  |  |  |
| Lava Ceram   | 3M Espe, Seefeld, Germany                  |  |  |  |  |  |
| Rondo Zirconia   | Nobel Biocare, Gothenborg, Sweden          |  |  |  |  |  |
| Triceram   | Dentaurum, Ispringen, Germany              |  |  |  |  |  |
| Vintage ZR   | Shofu, Kyoto, Japan                        |  |  |  |  |  |
| Vita VM9   | Vita, Bad Säckingen, Germany               |  |  |  |  |  |
| Zirox  | Wieland, Pforzheim, Germany                |  |  |  |  |  |
| Reflex   | Wieland, Pforzheim, Germany                |  |  |  |  |  |
| IPS d.sign   | Ivoclar-Vivadent, Schaan, Liechtenstein    |  |  |  |  |  |
| Vita VM13  | Vita, Bad Säckingen, Germany               |  |  |  |  |  |
| Veneering ceramics for th                              | e metal-ceramic technique are highlighted. |  |  |  |  |  |

# Table 1 Venering commiss used in the investigation

followed by a post hoc Scheffé test (SPSS Inc., Chicago, IL, USA; p < 0.05).

#### 2.1. Three-point flexural strength

Specimens with a final size of  $4 \pm 0.25$  mm in width,  $1.2 \pm 0.2$  mm in thickness and a length of at least 20 mm were produced.

The sample holder had a span between the two bearers of 15 mm. Supports and loading piston were steel knife edges, rounded to a radius of 0.8 mm. Load was applied at the midpoint of the specimens. The flexural strength was calculated according to the equation

$$\sigma = \frac{3 \text{Fl}}{2 b h^2}$$

where  $\sigma$  is the maximum center tensile stress (MPa), F the load at fracture (N), l the distance of the two supports (mm), b the width of the specimen (mm) and h is the height of the specimen (mm).

#### 2.2. Four-point flexural strength

Specimens with a final size of  $2.5\pm0.25\,mm$  in width,  $2.0\pm0.2\,mm$  in thickness and a length of at least 25 mm were used.

The sample holder had a span between the two bearers of 20 mm. The distance between the two loading pistons was 10 mm. Supports and both loading pistons were steel knife edges, rounded to a radius of 1.25 mm. The flexural strength was calculated according to the equation

$$\sigma = \frac{3Fd}{2bh^2}$$

where  $\sigma$  is the maximum center tensile stress (MPa), F the load at fracture (N), d the difference in the distance of the two supports and the distance of the two loading pistons (mm), b the width of the specimen (mm) and h is the height of the specimen (mm).

#### 2.3. Biaxial flexural strength

Disk-shaped specimens,  $12 \pm 0.2$  mm in diameter and  $1.2 \pm 0.2$  mm in height were prepared. The specimens were tested in a biaxial flexure jig with a piston on three balls design as described in the standard. The balls had a diameter of 3.2 mm and were arranged in an angle of  $120^{\circ}$  to each other on a circle of 10 mm in diameter. Loading at 1 mm/min was applied in the center of the specimen with a 1.5 mm diameter steel rod. Calculation of the biaxial flexural strength was performed with the following equation:

$$\sigma = \frac{-0.2387 F(X - Y)}{d^2}$$

where  $\sigma$  is the maximum center tensile stress (MPa), F the load at fracture (N),  $X = (1 + v) \ln(r_2/r_3)^2 + [(1 - v)/2] (r_2/r_3)^2$  and  $Y = (1 + v)[\ln(r_1/r_3)^2] + (1 - v)(r_2/r_3)^2$ .

In which, v is the Poisson's ratio,  $r_1$  the radius of the support circle (mm),  $r_2$  the radius of the loaded area (mm),  $r_3$  the radius of the specimen (mm) and d is the specimens thickness at the fracture origin (mm).

Poisson's ratio was taken as 0.25 for all ceramics according to the recommendation in the standard.

#### 3. Results

Means and respective standard deviations for three-point flexural strength, four-point flexural strength and biaxial flexural strength are shown in Table 3 and Fig. 1. For every ceramic the values of the three-point flexural strength were significantly higher than those of the four-point flexural

| Veneering<br>Ceramic | Pre Drying          |               | Heating          | Firing              | Holding       |
|----------------------|---------------------|---------------|------------------|---------------------|---------------|
|                      | Temperature<br>(°C) | Time<br>(min) | Rate<br>(°C/min) | Temperature<br>(°C) | Time<br>(min) |
| CerabienZR           | 600                 | 5             | 45               | 930                 | 1             |
| Creation ZI          | 450                 | 6             | 45               | 810                 | 1             |
| IPS e.max            | 400                 | 4             | 50               | 750                 | 1             |
| Initial ZR           | 400                 | 6             | 45               | 780                 | 1             |
| LavaCeram            | 450                 | 6             | 45               | 800                 | 1             |
| Rondo Zirconia       | 575                 | 5             | 45               | 925                 | 1             |
| Triceram             | 500                 | 6             | 55               | 760                 | 2             |
| Vintage ZR           | 650                 | 6             | 45               | 920                 | 1             |
| VM9                  | 500                 | 6             | 55               | 910                 | 1             |
| Zirox                | 575                 | 3             | 45               | 900                 | 2             |
| Reflex               | 575                 | 7             | 75               | 900                 | 2             |
| IPS d.sign           | 403                 | 6             | 60               | 869                 | 1             |
| VM13                 | 500                 | 6             | 55               | 880                 | 1             |

Vacuum was used until the final temperature was reached. Veneering ceramics for the metal-ceramic technique are highlighted.

| Veneering<br>Ceramic | 3-Point Flexural Strength<br>(MPa) |   | Biaxial Flexural Strength<br>(MPa) |     | 4-Point Flexural Strength<br>(MPa) |       |  |
|----------------------|------------------------------------|---|------------------------------------|-----|------------------------------------|-------|--|
| CerabienZR           | $77.8\pm8.7$                       | a | $69.1 \pm  4.8$                    | с   | $59.5\pm6.2$                       | g     |  |
| Vintage ZR           | $84.9 \pm 11.2$                    | a | $71.3 \pm 8.4$                     | с   | $64.8\pm6.3$                       | gh    |  |
| IPS e.max            | $85.7\pm20.5$                      | a | $73.2\pm10.4$                      | cd  | $69.2 \pm 5.1$                     | ghj   |  |
| Zirox                | $102.9\pm14.7$                     | a | $95.1 \pm 7.6$                     | e   | $71.7 \pm 4.4$                     | ghjk  |  |
| LavaCeram            | $90.0\pm9.0$                       | a | $86.1 \pm 7.0$                     | cde | $74.0 \pm  5.9$                    | ghjkl |  |
| Initial ZR           | $102.8\pm10.2$                     | a | $87.0 \pm  8.5$                    | cde | $80.2\pm6.3$                       | hjkl  |  |
| Creation ZI          | $98.7 \pm 17.2$                    | a | $92.3 \pm 9.6$                     | de  | $82.1 \pm 9.1$                     | jkl   |  |
| Rondo Zirconia       | $99.8 \pm 14.7$                    | a | $97.4 \pm 14.2$                    | e   | $82.6\pm10.1$                      | jkl   |  |
| Reflex               | $100.5\pm10.2$                     | a | $93.4\pm10.0$                      | de  | $82.7 \pm 8.5$                     | jkl   |  |
| IPS d.sign           | $93.3\pm13.5$                      | a | $95.5\pm7.8$                       | e   | $83.1 \pm 5.4$                     | jkl   |  |
| Triceram             | $105.5\pm12.4$                     | a | $101.4\pm10.5$                     | e   | $86.8\pm13.4$                      | kl    |  |
| VM9                  | $106.6\pm12.5$                     | a | $98.9 \pm 13.0$                    | e   | $89.2 \pm 9.5$                     | 1     |  |
| VM13                 | $149.4\pm20.5$                     | b | $141.2\pm11.6$                     | f   | $116.9\pm9.8$                      | m     |  |

Table 3 – Flexural strength values of the veneering ceramics (mean  $\pm$  S.D.), arranged in ascending order of the values for the four-point flexural strength

Identical letters following the values indicate homogeneous groups. Veneering ceramics for the metal-ceramic technique are highlighted.

strength. Statistical significant differences were found between three-point flexural strength and biaxial flexural strength for the following ceramics: Cerabien ZR, Initial ZR and Vintage ZR, while significant differences between biaxial flexural strength and four-point flexural strength occurred with Cerabien ZR, Lava Ceram, Rondo Zirconia, Triceram, Zirox and VM13. In Table 3, the homogeneous groups with no statistically significant differences between the different ceramics are marked. In the three-point flexure test, the strength values of the veneering ceramics for zirconia showed no statistically significant difference (group a). In the biaxial flexure test, three different homogeneous groups (c–e) of veneering ceramics for zirconia can be distinguished and in the four-point flexure test there were found five different groups (g, h, j, k, l) by statistical analysis. In the three-point flexure test, the values of the veneering ceramics for zirconia were similar to those of Reflex and IPS d.sign. In the biaxial flexure test, the flexure strengths of Cerabien ZR and Vintage ZR and in the four-point flexure test the flexure strengths of Cerabien ZR, Vintage ZR, IPS e.max, Zirox, Lava Ceram and Initial ZR were significantly lower than those of the veneering ceramics for the metal-ceramic technique. The flexural strength of VM13 in every case significantly exceeded those of the other ceramics investigated.



Fig. 1 - Flexural strength values and standard deviations of veneering ceramics.

Linear regression analysis revealed the following coefficients of determination:

three-point/four-point:  $R^2 = 0.89$ ,  $\sigma_{3-pt} = 1.24\sigma_{4-pt}$ ; three-point/biaxial:  $R^2 = 0.90$ ,  $\sigma_{3-pt} = 1.07\sigma_{biax}$ ; biaxial/four-point:  $R^2 = 0.92$ ,  $\sigma_{biax} = 1.16\sigma_{4-pt}$ .

### 4. Discussion

The results of this study revealed that the three-point flexural strength values of veneering ceramics for zirconia are in the same range as those of veneering ceramics for metal-ceramic systems. The regression analysis showed that the results of all three test methods are correlated. However, the three-point flexure test yielded the highest values. Compared to the fourpoint flexure test this difference was significant for all materials, compared to the biaxial flexure test only for three out of 13 ceramics. The biaxial flexure test in turn showed significantly higher values compared to the four-point test results for six out of 13 ceramics. But in general, it can be concluded that all three test designs provided the same relative order of the results. The four-point flexure test provided highest discrimination between the different ceramic materials, resulting in statistically significant differences between some veneering ceramics for zirconia and the control.

Similar biaxial flexural strength results as obtained in the present investigation are reported for leucite reinforced ceramics.<sup>21–23</sup> IPS d.sign showed a biaxial flexural strength of  $98.19 \pm 5.71 \text{ MPa}$ ,<sup>23</sup> which is comparable to the value measured in the present investigation (95.5  $\pm$  7.8 MPa). A further investigation employed biaxial flexure test and fourpoint flexure test.<sup>24</sup> Comparably low values for a body and an opaque ceramic for the metal-ceramic technique were found, but the relation between the results of both test methods was the same as in the present study. In another investigation, it is reported that IPS d.sign had a flexural strength in the threepoint, four-point and biaxial flexure strength test of  $124.3 \pm 12.4, 77.9 \pm 7.9$  and  $114.3 \pm 13.3$  MPa, respectively.<sup>20</sup> These values are quite high compared to the present investigation. Nevertheless, the authors also found a correlation between the three test methods, which was in the same order as in the present study. In a further study, it is reported that the three-point flexure strength of alumina was higher than that obtained in a biaxial flexure strength while this value was higher than the results obtained in a four-point flexure test, which again is in accordance with the present findings.<sup>25</sup>

The difference in the results of the three different test designs may be explained as follows. Flexural strength obtained with the four-point flexure test is generally lower because the probability to have a surface crack between the two loading pistons is higher than in the more limited area beneath the loading piston of a three-point flexure test. In the biaxial flexure test, the force is applied in the center of the specimen. Defects at the edges, which most probably lead to an early failure, are less effective. Nevertheless, the probability of a crack in the vicinity of the loading piston is higher than in the three-point flexure test because the loaded area is larger.<sup>19</sup> Consistent with Ban and Anusavice,<sup>24</sup> it can be concluded that for screening tests, for instance,

during the development of ceramics, the biaxial flexure test is most appropriate because preparation of the samples is easy, compared to the three- and four-point flexure tests. But, according to the present results, when a scientific approach is intended, the four-point flexure test should be preferred.

The fact that the strength of veneering ceramics for zirconia is in the same order as that of veneering materials for metal-ceramics could be interpreted in the sense that the strength of the veneering ceramics are not the limiting factor for the clinical long-term success of zirconia restorations. Nevertheless, compared to metal-ceramics excessive chipping is observed in clinical studies with zirconia restorations.<sup>3-4</sup> <sup>5</sup> To explain this effect, two aspects have to be considered. One aspect is the stress, built during cooling after firing of the veneering ceramic. In metal-ceramic systems, this stress may be at least partially relaxed by an elastic or plastic deformation of the substructure.<sup>15</sup> Especially, high-gold alloys show a low sag-resistance.<sup>16</sup> A zirconia substructure in contrast is rigid, which leads to higher stress formation. Hence, compared to metal-ceramics a higher flexural strength of the veneering ceramic is favorable to provide a high reliability of the veneer. The present investigation has shown that, depending on the test method and the brand, the flexural strength of veneering ceramics for zirconia is rather similar or even lower than that of veneering ceramics for the metal-ceramic technique. Therefore, the effort to improve the veneering ceramics for zirconia should be directed to the optimal adjustment of the thermal expansion and the increase of mechanical strength, which is in accordance with the appraisal of other authors.<sup>26</sup> A second point is the fact that in the oral cavity water exposure may cause hydrolysis of the Si-O-Si bonds, thus affecting the mechanical properties of the ceramic. Flexural strength values are obtained at ambient laboratory conditions. The increased failure rate of veneering ceramics for zirconia under humid conditions in the oral cavity may be attributed to a different chemical composition compared to ceramics for the metalceramic technique, resulting in a higher susceptibility for hydrolytic attack. Further investigations are scheduled to test this hypothesis.

# 5. Conclusion

Within the limitations of this *in vitro* study, the following conclusions can be drawn:

- (1) Four-point flexural strength values of all materials tested were significantly lower than those obtained with the three-point flexure test. The biaxial flexural strength in general ranged between the four-point flexural strength and the three-point flexural strength.
- (2) Strength values of veneering ceramics for zirconia are similar to those of veneering ceramics for the metalceramic technique.

# Acknowledgement

The materials were kindly provided by the respective manufacturers.

#### REFERENCES

- 1. Sturzenegger B, Fehér A, Lüthy H, Schärer P, Gauckler LJ. Reliability and strength of all-ceramic dental restorations fabricated by direct ceramic machining (DCM). International Journal of Computerized Dentistry 2001;4:89–106.
- Lüthy H, Filser F, Loeffel O, Schuhmacher M, Gauckler LJ, Hämmerle CHF. Strength and reliability of four unit allceramic posterior bridge. *Dental Materials* 2005;21:930–7.
- Vult von Steyern P, Carlson P, Nilner K. All-ceramic fixed partial dentures designed according to the DC-Zirkon technique. A 2-year clinical study. *Journal of Oral Rehabilitation* 2005;32:180–7.
- Sailer I, Fehér A, Filser F, Lüthy H, Gauckler LJ, Schärer P, et al. Prospective clinical study of zirconia posterior fixed partial dentures: 3-year follow-up. Quintessence International 2006;37:685–93.
- Sailer I, Fehér A, Filser F, Gauckler LJ, Lüthy H, Hämmerle CHF. Five-year clinical results of zirconia frameworks for posterior fixed partial dentures. *International Journal of Prosthodontics* 2007;20:383–8.
- al-Shehri SA, Mohammed H, Wilson CA. Influence of lamination on the flexure strength of dental castable ceramic. *Journal of Prosthetic Dentistry* 1996;76:23–8.
- Isgro G, Pallav P, van der Zel JM, Feilzer AJ. The influence of the veneering porcelain and different surface treatments on the biaxial flexure strength of a heat-pressed ceramic. *Journal of Prosthetic Dentistry* 2003;90:465–73.
- 8. De Jager N, Pallav P, Feilzer AJ. The influence of design parameters on the FEA-determined stress distribution in CAD-CAM produced all-ceramic crowns. *Dental Materials* 2005;**21**:242–51.
- Aboushelib MN, de Jager N, Kleverlaan CJ, Feilzer AJ. Microtensile bond strength of different components of core veneered all-ceramic restorations. *Dental Materials* 2005;21:984–91.
- Drummond JL, King TJ, Bapna MS, Koperski RD. Mechanical property evaluation of pressable restorative ceramics. *Dental* Materials 2000;16:226–33.
- Luthardt RG, Sandkuhl O, Reitz B. Zirconia-TZP and alumina—advanced technologies for the manufacturing of single crowns. The European Journal of Prosthodontics and Restorative Dentistry 1999;7:113–9.
- Aboushelib MN, de Jager N, Kleverlaan CJ, Feilzer AJ. Microtensile bond strength of different components of core veneered all-ceramic restorations. Part II: Zirconia veneering ceramics. Dental Materials 2006;22:857–63.
- 13. Al-Dohan HM, Yaman P, Dennison JB, Razzoog ME, Lang BR. Shear strength of core-veneer interface in

bi-layered ceramics. Journal of Prosthetic Dentistry 2004;**91**:349–55.

- Bagby M, Marshall SJ, Marshall GW. Metal ceramic compatibility: a review of the literature. *Journal of Prosthetic Dentistry* 1990;63:21–5.
- 15. Anusavice KJ, Carroll JE. Effect of incompatibility stress on the fit of metal-ceramic crowns. *Journal of Dental Research* 1987;**66**:1341–5.
- Fischer J, Fleetwood PW, Baltzer N. Thermal creep analysis of precious metal alloys for the ceramic-fused-to-metal technique. Journal of Biomedical Material Research (Applied Biomaterials) 1999;48:258–64.
- 17. Pjetursson BE, Sailer I, Zwahlen M, Hämmerle CHF. A systematic review of the survival and complication rates of all-ceramic and metal-ceramic reconstructions after an observation period of at least 3 years. Part I: single crowns. Clinical Oral Implant Research 2007;18:73–85.
- Pjetursson BE, Sailer I, Zwahlen M, Hämmerle CHF. A systematic review of the survival and complication rates of all-ceramic and metal-ceramic reconstructions after an observation period of at least 3 years. Part II: fixed partial dentures. Clinical Oral Implant Research 2007; 18:86–96.
- Anusavice KJ, Kakar K, Ferree N. Which mechanical and physical testing methods are relevant for predicting the clinical performance of ceramic-based dental prostheses? *Clinical Oral Implants Research* 2007;18:218–31.
- Jin J, Takahashi H, Iwasaki N. Effect of test method on flexural strength of recent dental ceramics. Dental Materials Journal 2004;23:490–6.
- Shareef MY, Van Noort R, Messer PF, Piddock V. The effect of microstructural features on the biaxial flexural strength of leucite reinforced glass ceramics. *Journal of Materials Science* Materials in Medicine 1994;5:113–8.
- Cattell MJ, Clarke RL, Lynch EJR. The biaxial flexural strength and reliability of four dental ceramics—part II. Journal of Dentistry 1997;25:409–14.
- Sinmazisik C, Övecoglu ML. Physical properties and microstructural characterization of dental porcelains mixed with distilled water and modeling liquid. *Dental Materials* 2006;22:735–45.
- Ban S, Anusavice KJ. Influence of test method on failure stress of brittle dental materials. *Journal of Dental Research* 1990;69:1791–9.
- Shetty DK, Rosenfield AR, Duckworth WH, Held PR. A biaxial-flexure test for evaluating ceramic strengths. *Journal* of the American Ceramic Society 1983;66:36–42.
- Aboushelib MN, de Jager N, Kleverlaan CJ, Feilzer AJ. Effect of loading method on the fracture mechanics of two layered all-ceramic restorative systems. *Dental Materials* 2007;23:952–9.

Abstract:

Objectives: The flexural strengths of veneering ceramics for zirconia were compared.Methods: With 10 different veneering ceramics for zirconia (test group) and three different veneering ceramics for the metal-ceramic technique (control group) three-point flexural strength and biaxial flexural strength according to ISO 6872: 1995 as well as four-point flexural strength according to EN 843-1: 2005 were measured (n = 10). Statistical analysis was performed with one-way ANOVA and post hoc Scheffe' test (SPSS, p < 0.05).

Results: For the test group, three-point flexural strength ranged between 77.8 +/- 8.7 and 106.6 +/- 12.5 MPa without any statistically significant differences, biaxial flexural strength between 69.1 +/- 4.8 and 101.4 +/- 10.5 MPa with three homogeneous groups and four-point

flexural strength between 59.5 +/- 6.2 and 89.2 +/- 9.5 MPa with five homogeneous groups. The control group showed three-point flexural strength values ranging from 93.3 +/- 13.5 to 149.4 +/- 20.5 MPa, biaxial flexural strength values from 93.4 +/- 10.0 to 141.2 +/- 11.6 MPa,

and four-point flexural strength values from 82.7 +/- 8.5 to 116.9 +/- 9.8 MPa. In every case, the results of the four-point flexure test were significantly lower than those obtained in the three-point flexure test. The three-point flexural strengths of the test group are similar to those of two ceramics of the control group. The flexural strength of one ceramic of the control group significantly exceeded the strengths of all other ceramics investigated. Conclusion: Three-point flexural strength values of veneering ceramics for zirconia are similar to those of veneering ceramics for the metal-ceramic technique. The four-point flexure test among all three tests showed highest discrimination between the different ceramic materials.