


B&D Zirconia, Kindest to Opposing Dentition



An adult molar with 4 distinct cusp tips was chosen for consistency in the density of a testing sample.



Then the tooth was sectioned into 4 equal parts.





Each different brand of zirconia sample was prepared according to Instructions for Use.




The sectioned tooth was fixed into a composite material in a jig to oppose each different brand of zirconia for comparative wear test.

B&D Zirconia, Kindest to Opposing Dentition


1.25 Million Cycles Simulates 5 Years of Chewing




B&D Dental
(Made in U.S.A.)
BeyondPlus™ Multi (A2)



Prominent Competitor 1
(Imported from Asia)
Multilayer zirconia (A2)



Prominent Competitor 2
(Made in U.S.A.)
Multilayer zirconia (A2)



Prominent Competitor 3
(Imported from Europe)
Multilayer zirconia (A2)









B&D Zirconia, Kindest to Opposing Dentition

B&D Dental
(Made in U.S.A.)
BeyondPlus™ Multi (A2)

Prominent Competitor 1
(Imported from Asia)
Multilayer zirconia (A2)

Prominent Competitor 2
(Made in U.S.A.)
Multilayer zirconia (A2)

Prominent Competitor 3
(Imported from Europe)
Multilayer zirconia (A2)



B&D Zirconia, Kindest to Opposing Dentition

1. ORIGIN Beyond Plus Multilayer (A2)



2. Prominent Competitor 1 (imported from Asia) Multilayer (A2)



3. Prominent Competitor 2 (Made in U.S.A.) Multilayer (A2)



4. Prominent Competitor 3 (imported from Europe) Multilayer (A2)

