

The True Next Generation High Performance Zirconia

- Anterior Translucency (>48%)
- Maximum Strength (>1250 MPa*)
- Eliminates Graying on Lighter Shades



3D-Master Shades & 16 Classic Shades

CHROMA™



The Perfect Culmination of Translucency and Strength, Beyond Plus™ zirconia is in a class of its own.

B&D Dental Technologies offers a uniquely universal zirconia that delivers on the promise of beauty and strength. **Beyond Plus™** is more than a true replacement for lithium disilicate such as ****e.max®**. It is stronger with equal or better aesthetics.

- Dentists no longer have to compromise strength for beauty. Their choice has been simplified.
- Labs benefit from a dense material that allows them to reproduce a sharp margin while still milling the 98mm zirconia disc easily (due to deliberate pre-sintering parameters).

In addition to the ultimate blend of translucency and strength, **Beyond Plus™** was developed to work with **CHROMA™** color liquid v.6.1 to provide consistent shading for optimal aesthetics.


*1100 MPa : up to 14 mm, 1250 MPa: over 16 mm **Not a trademark of B&D Dental



www.BnDdental.com
sales@BnDdental.com



800.255.2839



Translucency Comparison of **High Strength Zirconia (>1100 MPa*)**

"Beyond Plus™ is a **GAME CHANGER** for the industry. You get the **HT AESTHETICS** with the **STRENGTH** for multiple unit bridges. Finally we can make anterior as well as posterior bridges, with the strength needed."


Nadine Payne - Payne Oral Designs

*1100 MPa : up to 14 mm, 1250 MPa: over 16 mm

Translucency Comparison of **High Strength Zirconia**

"Beyond Plus™ is my **ZIRCONIA OF CHOICE** for **HIGHLY AESTHETIC CASES** that also demand **HIGH MPA STRENGTH**. It has great milling qualities and requires minimal staining after sintering.

Kevin Morgan, Morgan Dental Lab



** Up to 3 unit bridges. Larger restorations will require longer sintering cycles.

Sintering Schedule of **High Translucent, High Strength Zirconia (>1100 MPa*)**

"Beyond Plus™ simplifies my production with **ONE FAST** sintering cycle for **ALL** zirconia restorations. We have **CONSISTENT SUCCESS** and every shade is so close that it requires little or no stain."

Kirk Redding, Creative Crown Dental Lab

Comparison of **Strength & Translucency**

"Beyond Plus™ is **THE BEST ZIRCONIA ON THE MARKET**. It looks as good as any anterior zirconia, better than *e.max®, and is **STRONG** enough for large posterior bridges."

Jim Burke - JC'S Dental Designs LLC

